


PCR

TROUBLESHOOTING

THE ESSENTIAL GUIDE



Michael L. Altshuler



Caister Academic Press

TABLE OF CONTENTS

PREFACE	9
INTRODUCTION	11
APPEARANCES	13
• Unsatisfactory results of PCR	13
CAUSES AND ACTIONS	14
• Do not confront the problem	14
• The nature of PCR and its pathology	15
• Inadequate concentrations of ingredients	18
1. The template	18
2. Inadequate deoxynucleotidetriphosphates.....	20
3. Inadequate primers	20
4. Inadequate Taq	21
5. Inadequate Mg^{2+}	21
6. Suboptimal KCl concentration in the PCR buffer or the whole of the PCR buffer.....	21
• Inadequate quality of ingredients	23
1. DNA template.....	23
➤ Conversion of a DNA solution into a solid body	23
➤ PCR inhibitors.....	23
➤ Degraded template	23
➤ Verification of the purity of the template DNA by optical means.....	24
2. Poor water	24
3. Deoxynucleotidetriphosphates	24

4. Poor primers	25
➤ Primers may not be good in practice even if they are good in design	25
➤ For primer selection use only dedicated software.....	25
➤ Difference in T_m balanced by different primer concentrations.....	26
➤ Inefficient priming	26
5. Inadequate $MgCl_2$	26
6. Poor buffer	27
7. Poor Taq	27
8. The presence of PCR inhibitors.....	27
9. Substances that <u>do not</u> inhibit PCR	28
 • Inadequate storage of ingredients for the PCR reaction ..	29
1. The treacherous refrigerators.....	29
2. Templates, dNTPs, primers	29
3. Taq polymerase	31
4. $MgCl_2$	31
5. Buffer	31
 • The thermocycler	32
1. Conductivity of heat puts a limit to the mix composition	32
2. The ramp	33
3. Time and temperature.....	34
4. Dusty, greasy, or fluffy tube wells.....	36
5. Rapid evaporation.....	36
6. Suboptimal performance of the thermocycler or its particular wells	36
7. The tubes are poorly pressed down into the wells or they have got deformed.....	37

• Faulty target selection	38
• Incomplete DNA denaturation and dispersal	40
1. Template DNA	40
2. PCR fragments	41
3. Hairpins	42
• The Taq enzyme	43
1. Hurdles for Taq polymerase	43
➤ Stable hairpins in the template strand	43
➤ AT-rich areas	43
➤ GC-rich areas	43
➤ Alternating GC/AT-rich regions	43
2. Hot start	44
➤ The improved hot start.....	44
➤ The role of the reaction volume in the quasi hot start.....	45
3. Nonspecific binding of Taq to DNA	46
• Incomplete primer elongation or premature termination of DNA synthesis	48
1. Under-elongation of primers in the late PCR	48
2. Premature termination of the DNA synthesis.....	50
• Cosolvents or additives or enhancers	51
1. Helix-destabilizers.....	51
2. Helix-stabilizers.....	52
3. Substances that neutralise the PCR inhibitors	53
4. PCR enhancers with poorly understood mechanism of action.....	53
• Approaching the limit of the PCR sensitivity	54

• Agarose gel electrophoresis	57
1. The band diffuses and disappears	57
2. Short fragments of uneven length migrate down the gel without any separation	57
3. The band is invisible.....	57
4. The significance and the insignificance of the salt concentration in the compared samples.....	57
5. Bands smear due to their fast movement or the DNA overload	58
6. Dirty gel support.....	58
7. Dried well	58
8. DNA sticking in the gel well caused by inappropriate gel density	58
• Causes for specific nonspecifics and the false contamination	59
1. Chimeras	59
2. Allele dropout.....	59
3. Heteroduplexes	60
4. Primer multimers	61
5. Low resolution electrophoresis resulting in imprecise idea of the correct band position.....	61
6. Coincidence or the devil.....	61
• Mineral oil and wax	62
1. Mineral oil	62
2. Wax, paraffin or vaseline.....	62
• Primer-dimers and primer multimers	63
• Short PCR fragments versus long PCR fragments	67
• Avoiding accidents	68

CONCLUSIONS	69
• A few words to the novice	69
• A few words to a PCR adept	69
GLOSSARY	71
INDEX	77